

ESTHETIC CASE INVENTORY, DETAILED

1. Detailed laboratory instructions and narrative, to include:
 - Teeth treated.
 - Type of restorations (if established).
 - Reason for treatment.
 - Patient's goals and desires.
 - Relevant treatment information ie: future Tx, recent bleaching, health-compromised patient, etc..).
 - Inventory of items sent.
 - Two shades desired (taken before procedures) and existing/underlying.
 - Shade mapping as a complement to photos.
2. Photographic records for smile design: (minimum AACD standard views required)
 - View 1: Full face, frontal, 1:10.
 - View 2: Full smile, non-retracted frontal, 1:2.
 - View 5: Retracted frontal view, 1:2.
 - View 8: Retracted frontal view incisors (or treated teeth), 1:1.
 - Any other view that depicts relevant perspective.
3. Photographic records for shade:
 - High-quality slides, prints or high-resolution digital prints showing teeth to be matched and the closest shade guide together, at same distance from camera.
 - Must be taken at very beginning of appointment to avoid dehydration lightening.
 - The number on the shade tabs must be visible.
 - Take at least 3 photos at different angles to avoid reflections.
 - Opposing arch: same as above, depicting teeth of opposing arch.
 - Views of prepped teeth with closest stump shade.
4. Impression of the prepared arch:
 - Polyvinyl or polyether impression, so it may be poured several times.
 - Full arch.
 - All margins clearly visible.
 - All occlusal surfaces of posteriors clear, for correct mounting.
5. Opposing model:
 - Preferably taken in a resilient stable material, same as treated arch.
 - Full arch.
 - If alginate, then to be poured with vacuum-mixed stone.

6. Approved provisional model, patient-approved for form and function:
 - Full arch so it can be mounted in centric.
 - The most accurate and clear possible.
 - Must clearly capture gingival contours, embrasures.
 - Posterior occlusal surfaces clear and undistorted.
7. Pre-operative model:
 - Useful for diagnosis and reference.
8. Two types of bite records:
 - Between Preparations and Opposing
 - Centric bite, not wax or Duralay, but elastomeric material.
 - Ensure patient is not retruded or underclosed
 - Between Provisionals and Opposing (VERY IMPORTANT)
 - So study model is mounted correctly and followed.
9. Horizontal reference:
 - Midliner® with bubble level glued in place
 - (WE HAVE MIDLINERS AVAILABLE - PLEASE ASK US).
 - Stick bite (use plastic brush handle instead of Q-tip wood applicator).
10. Any other relevant item or record, ie:
 - Reference from a smile guide.
 - Cut-out from a magazine.
 - Patient's previous photograph.
 - X-rays, etc...

